



- d. Remove the sample from the oven and carefully transfer it to an agate mortar. Grind the sample to pass a standard 180- $\mu$ m sieve, and return it to the vacuum oven at 100 kPa for approximately 1 h.
2. Solid Sample Preparation:
    - a. Depending on their composition, solid samples may be handled in several ways. As an example, if the sample is a soluble resinous material, it may be possible to dissolve a portion in solvent and evaporate the solution directly on an alkali halide disk.
    - b. Non-soluble solids should be ground into pieces and dried overnight in an aluminum dish in the vacuum oven at 60°C and 100 kPa.
  3. Remove the sample from the oven. Grind 2 mg of the sample with 250 mg of infrared grade potassium bromide, KBr, until it forms a uniform mixture. This breaks up any lumps of KBr and provides a preliminary mix to the specimen.
  4. Transfer the mixed material, plus one steel ball, to a Wig-L-Bug capsule and vibrate the sample for 30 to 60 s.
  5. Being careful that the steel ball is not also transferred, transfer the powdered pellet specimen to the evacuable press, and follow the manufacturer's instructions to prepare a suitable disk. Thoroughly clean the die after each use. Be careful to avoid damaging the polished die faces.
  6. Place the disk in an infrared spectrophotometer and collect the transmittance spectrum.

Test results are used for comparison purposes only. Each spectrum is compared with samples run previously. Two materials are considered similar if all of

the absorption peaks match as to wave length and relative magnitude.

## PART 2. SAFETY AND HEALTH

This method may involve hazardous materials, operations, and equipment. This method does not purport to address all the safety problems associated with its use. It is the responsibility of the user of this method to consult and establish appropriate safety and health practices and determines the applicability of regulatory limitations prior to use.

Prior to handling, testing or disposing of any of waste materials, testers are required to read the Caltrans Laboratory Safety Manual. This manual contains information on general safety principles, standard operating procedures, protective apparel, disposal of materials and how to handle spills, accidents, emergencies, etc. Users of this method do so at their own risk.

### REFERENCES:

Caltrans Laboratory Safety Manual  
California Test 416

End of Test (California Test 416 contains 2 pages)